

Solenoid Relay Valve Kit

Part No. 59947-05K

- Valve
- Sleeve

NOTE: Also fits AF23/33 & RE5F22A units.



1. Disassembly

- a. Remove and save OE retaining clip for reuse.
- b. Remove and discard OE sleeve and valve.

2. Installation & Assembly

- a. Install Sonnax valve, followed by Sonnax sleeve.
- b. Reinstall OE retaining clip.

3. Final Testing

Vacuum testing at the port(s) indicated holds the recommended minimum 18 in-Hg.

